

Confirmation Flippers

20D

Over →

Optego's confirmation flippers have a molded plastic frame with chemically hardened glass lenses. All flippers come standard with a beige frame unless an alternative color is specified or requested. Made in the USA.

Standard Diagnostic Flipper Set

This set includes 6 flippers with colored frames and a flipper base.

Base included

+/- 0.25 D - purple

+/- 0.50 D - green


+/- 0.75 D - yellow

+/- 1.00 D - red

+/- 1.50 D - blue

+/- 2.00 D - orange

93450



OEP Chairside Diagnostic Flipper Set

This set includes 6 flippers with colored frames, with 3 plus/plus flippers and 3 of the most common plus/minus flippers.

Base included

+/- 0.25 D - red

+/- 0.50 D - orange

+0.50 / +0.75 D - yellow

+1.00 / +1.25 D - green

+1.50 / +1.75 D - blue

+2.00 / -2.00 D - purple

93451



Flipper Base

93459

Contact Lens Flippers

Side 1: +0.25 / -0.25 D

Side 2: +0.50 / -0.50 D

93455

Computer Home Vision Therapy Flippers

Flipper A

93456

Side 1: +0.75 / -1.50 D // Side 2: -2.50 / +1.25 D

Flipper B

93457

Side 1: +1.75 / -3.50 D // Side 2: -4.00 / +2.00 D

Flipper C

93458

Side 1: +2.25 / -4.50 D // Side 2: -5.00 / +2.50 D

Stock Spherical Flippers

Two pairs of spherical lenses: one pair of plus lenses (black labels) and one pair of minus lenses (red labels).

Dpt.	Prod. #	Dpt.	Prod. #
+/- 0.25 D	94200	+1.00 D/-2.00 D	94253
+/- 0.50 D	94201	+1.00 D/-3.00 D	94254
+/- 0.75 D	94202	+1.25 D/-2.50 D	94248
+/- 1.00 D	94203	+1.50 D/-3.00 D	94249
+/- 1.25 D	94204	+2.00 D/-2.50 D	94255
+/- 1.50 D	94205	+2.00 D/-3.00 D	94256
+/- 1.75 D	94206	+2.00 D/-3.50 D	94257
+/- 2.00 D	94207	+2.00 D/-4.00 D	94258
+/- 2.25 D	94208	+2.00 D/-4.50 D	94259
+/- 2.50 D	94209	+2.00 D/-5.00 D	94260
+/- 2.75 D	94210	+2.00 D/-6.00 D	94261
+/- 3.00 D	94211	+2.50 D/-3.00 D	94262
+/- 3.50 D	94212	+2.50 D/-4.00 D	94263
+0.50 D/+1.00 D	94250	+2.50 D/-4.50 D	94264
+0.50 D/-1.00 D	94251	+2.50 D/-5.00 D	94265
+0.75 D/+1.25 D	94252		

Stock Prism Flippers

Two pairs of prism lenses: one pair base "in" and one pair base "out". Diopter and base orientation (black labels).



Dpt.	Prod. #	Dpt.	Prod. #
1^in / 1^out	94270	2^in / 4^out	94285
2^in / 2^out	94271	3^in / 12^out	94286
3^in / 3^out	94272	4^in / 6^out	94287
4^in / 4^out	94273	4^in / 8^out	94288
5^in / 5^out	94274	4^in / 12^out	94289
6^in / 6^out	94275	5^in / 10^out	94290
7^in / 7^out	94276	6^in / 8^out	94291
8^in / 8^out	94277	6^in / 10^out	94292
9^in / 9^out	94278	6^in / 12^out	94293
10^in / 10^out	94279	8^in / 10^out	94294
12^in / 12^out	94280	8^in / 12^out	94295
14^in / 14^out	94281	8^in / 14^out	94296
16^in / 16^out	94282	8^in / 16^out	94297
18^in / 18^out	94283	10^in / 12^out	94298
20^in / 20^out	94284	10^in / 20^out	94299

Custom Spherical Flippers/ Prism Flippers (beige)
For any prism power combination that is not a stock power listed above (max. 20^ each lens), please visit optegousa.com, call 1-866-678-6723 or use our Custom Flipper Order Form on the back.

OPTEGO VISION USA
4525 South Coach Dr
Tucson, AZ 85714
Phone 1-866-678-6723
Fax 1-520-209-1586
email@optego.com
optegousa.com

OPTEGO VISION ASIA
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
optego.com.sg

Date		PO #	
Company Name		Contact Name	
Shipping Address			
City	State	ZIP	
Phone:		Email:	



Plus Sphere Lenses: +0.12 / +0.25 / +0.37 / +0.50 / +0.62 / +0.75 / +1.00 / +1.25 / +1.50 / +1.75 / +2.00 / +2.25 / +2.50 / +2.75 / +3.00 / +3.25 / +3.50 / +3.75 / +4.00 / +4.50 / +5.00

Minus Sphere Lenses: -0.12 / -0.25 / -0.37 / -0.50 / -0.62 / -0.75 / -1.00 / -1.25 / -1.50 / -1.75 / -2.00 / -2.25 / -2.50 / -2.75 / -3.00 / -3.25 / -3.50 / -3.75 / -4.00 / -4.50 / -5.00 / -5.50 / -6.00 / -6.50 / -7.00 / -8.00

Type	1st Lens	2nd Lens	3rd Lens	4th Lens	Handle	Qty
Sphere					☐ Right ☐ Left	
Sphere					☐ Right ☐ Left	
Sphere					☐ Right ☐ Left	
Sphere					☐ Right ☐ Left	
Sphere					☐ Right ☐ Left	

Prism Lenses: 0.25^Δ / 0.5^Δ / 1^Δ / 1.5^Δ / 2^Δ / 2.5^Δ / 3^Δ / 3.5^Δ / 4^Δ / 4.5^Δ / 5^Δ / 6^Δ / 7^Δ / 8^Δ / 9^Δ / 10^Δ / 12^Δ / 14^Δ / 16^Δ / 18^Δ / 20^Δ

Type	1st Lens	2nd Lens	3rd Lens	4th Lens	Handle	Qty
Prism					☐ Right ☐ Left	
Lens Base	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down		
Prism					☐ Right ☐ Left	
Lens Base	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down		
Prism					☐ Right ☐ Left	
Lens Base	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down		
Prism					☐ Right ☐ Left	
Lens Base	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down	☐ in ☐ out ☐ up ☐ down		