

Occluders / Maddox Rods

Handheld Monocular Occluders

Black Length: 9"

93115

White Length: 9"

93116



Handheld Double-ended Occluders

Black Occluder and Maddox Lens
Length: 9-5/8"

93141



Red Lens and Maddox Lens
Length: 9-5/8"

93147



Handheld Single Pinhole Occluder

9313

Black Long handle: 9"



Handheld Maddox Rod

Handheld Maddox Rod
Long handle
Length: 9"

9314



Handheld Multiple Pinhole Occluder

9312


Black Long handle: 9"



Spielmann Translucent Occluder

93501

Translucent Occluder
Length: 8-1/2"
(23.5 cm)



Lorgnette Pinhole Occluder

93101



With multiple pinhole attachment
17 holes (1.5 mm diam.)

The multiple pinhole attachment
will stay in its "upright" position
due to an improved design.

Strong, self-formed rivet for long
durability.

The occluder is designed to provide occlusion while at
the same time giving the practitioner the opportunity to
observe any movements of the occluded eye.

OPTEGO VISION USA
4525 South Coach Dr
Tucson, AZ 85714
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
optegousa.com

OPTEGO VISION ASIA
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
optego.com.sg