

Hypo-allergenic Occlusion Eye Patches

Ortopad is an occlusion patch used in orthoptic treatment of amblyopia and strabismus.

The eye patches are made from a soft, air-permeable, eco-friendly bamboo material with a light protection inlay and a special hypo-allergenic "hotmelt" adhesive. They are latex-free and contain no preservatives, which means they can be used on a daily basis, even with sensitive skin..

The shape of the occlusion patch adapts perfectly to the area around the eye. Ortopad® patches have two nasal slits for a better fit.

- 50 patches per box, individually wrapped
- Ortopad® Boys and Ortopad® Girls (Junior, Medium and Regular)
- All boxes with Ortopad® patches (except white and beige) have 5 different patch designs
- Each box of White and Beige Ortopad® patches contains 24 colored and 24 blank stickers. Fun for children to color and to decorate the patch
- One instruction sheet per box
- Ortopad patches are latex free
- Complimentary samples
- Conforms to European Standards (CE)
- Made in Italy



We make patching fun!!

ORTOPAD PATCHES 50/box

Junior (up to 2 years)	2-5/8" x 2" / 67 x 50 mm
Medium (2 - 4 years)	3" x 2-1/8" / 76 x 54 mm
Regular (over 4 years)	3-5/8" x 2-1/4" / 85 x 59 mm

ORTOPAD GIRLS



ORTOPAD BOYS



ORTOPAD Girls/Boys 50/box	Prod #
Girls Junior	33000013
Girls Medium	33000014
Girls Regular	33000015
Boys Junior	33000016
Boys Medium	33000017
Boys Regular	33000018

OPTEGO VISION USA
4525 South Coach Dr
Tucson, AZ 85714
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
optegousa.com

OPTEGO VISION ASIA
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
www.optego.com.sg

ORTOPAD SOFT GIRLS



Ortopad SOFT has extra soft textured accents

ORTOPAD SOFT BOYS



ORTOPAD SOFT Girls/Boys 50/box	Prod #
Girls Junior	33000019
Girls Medium	33000020
Girls Regular	33000021
Boys Junior	33000022
Boys Medium	33000023
Boys Regular	33000024

ORTOPAD FUN GIRLS



ORTOPAD FUN BOYS



ORTOPAD FUN Girls/Boys 50/box	Prod #
Girls Junior	33000993
Girls Medium	33000994
Girls Regular	33000995
Boys Junior	33000773
Boys Medium	33000774
Boys Regular	33000075

ORTOPAD FUN PACK



ORTOPAD FUN PACK 50/box	Prod #
Junior	33000010
Medium	33000011
Regular	33000012

Each box of White and Beige Ortopad® patches contains 24 colored and 24 blank stickers. Fun for children to color and to decorate the patch.

ORTOPAD WHITE



ORTOPAD BEIGE



ORTOPAD White/Beige 50/box	Prod #
White Junior	33000004
White Medium	33000005
White Regular	33000006
Beige Junior	33000007
Beige Regular	33000009