

Prism Blocks

Prism Blocks

Single universal, acrylic. Height: 37 mm. The dioptric power of each prism block has been laser engraved. Superior surface finish!



Dpt.	Prod. #	Dpt.	Prod. #
Prism 0.5 D	936700	Prism 14 D	936714
Prism 1 D	936701	Prism 15 D	936715
Prism 2 D	936702	Prism 16 D	936716
Prism 3 D	936703	Prism 18 D	936718
Prism 4 D	936704	Prism 20 D	936720
Prism 5 D	936705	Prism 25 D	936725
Prism 6 D	936706	Prism 30 D	936730
Prism 7 D	936707	Prism 35 D	936735
Prism 8 D	936708	Prism 40 D	936740
Prism 9 D	936709	Prism 45 D	936745
Prism 10 D	936710	Prism 50 D	936750
Prism 12 D	936712		

Prism Blocks on a handle ("Stick Prisms")

A very effective way of examining patients. Easier to hold in front of patient's eye. The handle is made of knurled aluminum (length: 4.1"/10.5 cm) and is permanently attached to the prism.



Dpt.	Prod. #	Dpt.	Prod. #
Prism 0.5 D	936900	Prism 14 D	936914
Prism 1 D	936901	Prism 15 D	936915
Prism 2 D	936902	Prism 16 D	936916
Prism 3 D	936903	Prism 18 D	936918
Prism 4 D	936904	Prism 20 D	936920
Prism 5 D	936905	Prism 25 D	936925
Prism 6 D	936906	Prism 30 D	936930
Prism 7 D	936907	Prism 35 D	936935
Prism 8 D	936908	Prism 40 D	936940
Prism 9 D	936909	Prism 45 D	936945
Prism 10 D	936910	Prism 50 D	936950
Prism 12 D	936912		

Prism Blocks Sets

Prism Blocks Sets

Prism Blocks Sets in a compact polymer carrying case. The box opens flat and fits in a drawer. Superior optical quality!

Size: 8-7/8" W x 8.0" D x 2.0" H (22.5 x 20 x 5 cm)

A set of 22 acrylic prisms in polymer case 92352
0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50 D and separate red lens.



A set of 16 acrylic prisms in polymer case 92353
0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45, 50 D



All prisms are manufactured following Optego's proprietary ClearPrism® Technology Protocol.

OPTEGO VISION USA
4525 South Coach Dr
Tucson, AZ 85714
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
optegousa.com

OPTEGO VISION ASIA
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
optego.com.sg