

Plastic Eye Patches


Plastic Eye Patches
3" x 2" with elastic band
for easy adjustment.

Rigid plastic
black 1/pkg 98570
black 10/pkg 98540

Flexible plastic
black 1/pkg 98577
black 10/pkg 98541


Flexible Plastic Shield with edge cushion

Flexible plastic shield (#98577) with non-allergenic silicone overmold edge cushion. The removable edge cushion can be thoroughly cleaned with alcohol. Very comfortable for patients during visual field examinations.



Available overmold colors

Blue	98560
Gray	98565
Green	98562
Purple	98564
Red	98566
Yellow	98568



Fox Aluminum Eye Shields

97880

Assorted colors, 10/pkg

Fox Aluminum Eye Shields with removable, malleable silicone garter. Silicone garter can be removed, cleaned and easily repositioned on shield.



Pressure Patch (Optego)

98510


Non-sterile 1/bag
Compare to PressPatch™

The pressure patch consists of a plastic shield with a foam pad, kept in place over the eye with an adjustable elastic strap.



For Drape Support, see Miscellaneous section

Eye Shields (vented)



Clear Eye Shields (vented)

Polycarbonate, clear 10/pkg
Polycarbonate, clear 50/box

97836
97839

Our polycarbonate shields are optically clear, without distortion. Universal, for left and right eye. Peripheral vents for patient comfort.

Opaque Eye Shields (vented)

Plastic, white 10/pkg
Universal, for left and right eye.

98579

Karickhoff Visual Field Fatigue Diffuser 98571

3/pkg

White plastic, with elastic band. Makes visual fields more accurate and comfortable by allowing the same amount of light to enter both eyes. Eliminates "bowl blackout" of the tested eye.



Ortolux® Eye Bubble Moisture Chamber

Small adhesive (95 x 63 mm)	20/pkg	94122
Large adhesive (111 x 79 mm)	20/pkg	94124
Large adhesive (111 x 79 mm)	1/pkg	94121

Clear dome for clear central vision with latex free adhesive edge. Helps to prevent the cornea from drying out. Ideal product for patients with Bell's Palsy.



Ortolux® Air Eye Shield™ (vented)

Small adhesive (16 x 24 mm)	20/box	94123
Large adhesive (22 x 35 mm)	20/box	94120

Chamber size: 56 x 37 mm
Universal, for left and right eye
Clear, polycarbonate vented eye shield with latex free adhesive edge; no need for separate tape.



OPTEGO / USA
3391 E Hemisphere Loop
Tucson, AZ 85706
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
www.optegousa.com

OPTEGO / CANADA
3600B Laird Rd, Unit 3
Mississauga ON L5L 6A7
Phone 1 416-780-1289
Toll Free 1 877-678-3464
Fax 1 877-471-6001
orders@optego.com
www.optego.com

OPTEGO / SINGAPORE
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
www.optego.com.sg