

Surgical Drape Support / Retractors 98501

Disposable Drape Support 20/pkg




The drape support keeps the drape away from the mouth and nose during surgical procedures: breathing becomes much easier. 10-5/8" x 1-5/8"

Barraquer Wire Speculums

Barraquer Wire Speculum
Open blades with bar, reusable

8201

Adult
Blade Width: 15 mm
Blade Spread: 17 mm



Barraquer Wire Speculum
Solid blades, reusable

8207

Adult (matte finish)
Blade Width: 14 mm
Blade Spread: 19 mm



Phorofter Near Vision Rod 93692



Holds Phorofter clip with near vision card

Optokinetic Tape

OKN Tape with black/white stripes 93390
OKN Tape with red/white stripes 9339



Striped fabric with black and white or red and white stripes.
Rollable and foldable.
Size : 28-5/8" long x 3-3/8" tall.
Stripe size: 2-5/16"

Alcohol Dispenser 9949

Liquid well pump
7.6 fl.oz. (225 ml)

Excellent for dispensing isopropyl alcohol. Bottle clearly marked with 'Alcohol' in red.
The clear bottle allows for fluid level visibility. Unique push design, preventing cross contamination.



OPTEGO / USA
3391 E Hemisphere Loop
Tucson, AZ 85706
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
www.optegousa.com

OPTEGO / CANADA
3600B Laird Rd, Unit 3
Mississauga ON L5L 6A7
Phone 1 416-780-1289
Toll Free 1 877-678-3464
Fax 1 877-471-6001
orders@optego.com
www.optego.com

OPTEGO / SINGAPORE
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
www.optego.com.sg