


Prism Blocks

All prisms are manufactured following Optego's proprietary ClearPrism® Technology Protocol.

Prism Blocks

Single universal, acrylic. Height: 37 mm. The dioptic power of each prism block has been laser engraved. Superior surface finish!



| Dpt. | Prod. # | Dpt. | Prod. # |
|-------------|---------|------------|---------|
| Prism 0.5 D | 936700 | Prism 14 D | 936714 |
| Prism 1 D | 936701 | Prism 15 D | 936715 |
| Prism 2 D | 936702 | Prism 16 D | 936716 |
| Prism 3 D | 936703 | Prism 18 D | 936718 |
| Prism 4 D | 936704 | Prism 20 D | 936720 |
| Prism 5 D | 936705 | Prism 25 D | 936725 |
| Prism 6 D | 936706 | Prism 30 D | 936730 |
| Prism 7 D | 936707 | Prism 35 D | 936735 |
| Prism 8 D | 936708 | Prism 40 D | 936740 |
| Prism 9 D | 936709 | Prism 45 D | 936745 |
| Prism 10 D | 936710 | Prism 50 D | 936750 |
| Prism 12 D | 936712 | | |

Prism Blocks on a handle ("Stick Prisms")

A very effective way of examining patients. Easier to hold in front of patient's eye. The handle is made of knurled aluminum (length: 4.1"/10.5 cm) and is permanently attached to the prism.




| Dpt. | Prod. # | Dpt. | Prod. # |
|-------------|---------|------------|---------|
| Prism 0.5 D | 936900 | Prism 14 D | 936914 |
| Prism 1 D | 936901 | Prism 15 D | 936915 |
| Prism 2 D | 936902 | Prism 16 D | 936916 |
| Prism 3 D | 936903 | Prism 18 D | 936918 |
| Prism 4 D | 936904 | Prism 20 D | 936920 |
| Prism 5 D | 936905 | Prism 25 D | 936925 |
| Prism 6 D | 936906 | Prism 30 D | 936930 |
| Prism 7 D | 936907 | Prism 35 D | 936935 |
| Prism 8 D | 936908 | Prism 40 D | 936940 |
| Prism 9 D | 936909 | Prism 45 D | 936945 |
| Prism 10 D | 936910 | Prism 50 D | 936950 |
| Prism 12 D | 936912 | | |

Prism Blocks Sets

Magnetic Prism Block Set (in polymer box)

92354

Each block prism contains a small magnet and 'snaps on' to a magnetic handle. Length: 4" (10 cm)
A set of 22 prisms. Two magnetic handles
0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50 D (no red lens)



Prism Blocks Sets

In a compact polymer case with carrying handle. The box opens flat and fits in a drawer. Superior optical quality!
Size: 8-7/8" W x 8.0" D x 2.0" H (22.5 x 20 x 5 cm)

A Set of 22 Acrylic Prisms in polymer box

92352

0.5, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 12, 14, 16, 18, 20, 25, 30, 35, 40, 45, 50 D and separate red lens.



A Set of 16 Acrylic Prisms in polymer box

92353

0.5, 1, 2, 4, 6, 8, 10, 12, 15, 20, 25, 30, 35, 40, 45, 50 D



OPTEGO / USA
3391 E Hemisphere Loop
Tucson, AZ 85706
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
www.optegousa.com

OPTEGO / CANADA
3600B Laird Rd, Unit 3
Mississauga ON L5L 6A7
Phone 1 416-780-1289
Toll Free 1 877-678-3464
Fax 1 877-471-6001
orders@optego.com
www.optego.com

OPTEGO / SINGAPORE
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
www.optego.com.sg