

## Retinoscopy Racks

93610

Two Paddle School Bus Retinoscopy Racks



Two paddle retinoscopy racks with glass lenses, shaped like a school bus.  
One paddle with plus lenses and a second paddle with minus lenses (+/- 0.5 D to +/- 15.0 D)  
Comes with protective case

## Retinoscopy Rack

93608

One Paddle ABCD Retinoscopy Rack



One paddle with glass lenses, shaped like a schoolbus.  
This paddle has plus lenses 1, 2, 3, 4, 5, 6, 7, 8, 9, 10  
one frosted lens and one -5 lens  
Inspired by Dr. Robert Arnold at Alaska Blind Child  
Discovery (ABCD)

OPTEGO / USA  
3391 E Hemisphere Loop  
Tucson, AZ 85706  
Phone 866-678-6723  
Fax 520-209-1586  
email@optego.com  
www.optegousa.com

OPTEGO / CANADA  
3600B Laird Rd, Unit 3  
Mississauga ON L5L 6A7  
Phone 1 416-780-1289  
Toll Free 1 877-678-3464  
Fax 1 877-471-6001  
orders@optego.com  
www.optego.com

OPTEGO / SINGAPORE  
60 Paya Lebar Rd #11-33  
Paya Lebar Square  
Singapore 409051  
Phone +65 6742 5940  
asia@optego.com  
www.optego.com.sg