

Miscellaneous Products

Surgical Drape Support 98501

Disposable Drape Support 20/pkg



The drape support keeps the drape away from the mouth and nose during surgical procedures: breathing becomes much easier. 10-5/8" x 1-5/8"

Phoropter Near Vision Rod 93692


detail of rod



Holds Phoropter clip with near vision card

Red-Green Glasses - reversible 9397

Red-Green reversible glasses
Chrome frame



Red-Green Glasses - adult 9374

1 pair
2 pairs

9374
93748

Plastic, 1 pair
adult size



Red lens on right side - green lens on left

Plastic, 2 pairs
adult size



1 pair red lens on right - green lens on left
1 pair red lens on left - green lens on right

Red-Green Glasses - pediatric size 9375

Plastic, 1 pair,
pediatric size
Frame color may vary



Barraquer Wire Speculum 8207

Barraquer Wire Speculum
Solid blades, reusable
Adult (matte finish)
Blade Width: 14 mm
Blade Spread: 19 mm



Alcohol Dispenser 9949

Liquid well pump
7.6 fl.oz.(225 ml)

Excellent for dispensing isopropyl alcohol. The clear bottle allows for fluid level visibility. Unique push design, preventing cross contamination.




Eye Model with Muscles 93901



Detailed 2-inch eye model showing 6 extra-ocular muscles, lacrimal gland and full range of parts of the eye. Front hinged. Fully assembled version mounted on a solid plastic base.

Eye Model Small 93900



Highly detailed eye model, dissectable into 6 parts: top section, vitreous body, bottom section, lens, iris and cornea. Mounted on a plastic stand. 3-1/2" x 3-1/2" x 4-3/8"

Eye Model Large 93902



Highly detailed eye model, dissectable into 6 parts: top section, vitreous body, bottom section, lens, iris and cornea. Mounted on a plastic stand. 5-1/2" x 5-1/2" x 6 1/4"

OPTEGO / USA
3391 E Hemisphere Loop
Tucson, AZ 85706
Phone 866-678-6723
Fax 520-209-1586
email@optego.com
www.optego.com

OPTEGO / CANADA
3600B Laird Rd, Unit 3
Mississauga ON L5L 6A7
Phone 1 416-780-1289
Toll Free 1 877-678-3464
Fax 1 877-471-6001
orders@optego.com
www.optego.com

OPTEGO / SINGAPORE
60 Paya Lebar Rd #11-33
Paya Lebar Square
Singapore 409051
Phone +65 6742 5940
asia@optego.com
www.optego.com.sg